

THE MOST IN-DEMAND  
AESTHETIC TREATMENTS IN

ONE

WORKSTATION



VENUS VERSA

# A REVOLUTIONARY

ADVANCE IN

**VERSATILITY**



## **PHOTOREJUVENATION**

---

Sun damage reduction,  
color correction, vascular lesion  
reduction, skin tone improvement

## **HAIR REMOVAL**

---

Treats all body parts,  
including underarms,  
face, back and bikini area

## **ACNE TREATMENT**

---

Active acne and  
inflammation reduction

## **SKIN RESURFACING**

---

Acne scar reduction  
pore minimization,  
deep wrinkle reduction,  
texture improvement

## **BODY CONTOURING**

---

Cellulite reduction and  
circumferential reduction

## **ANTI-AGING**

---

Skin tightening and  
wrinkle reduction

# IPL APPLICATORS & OUR TECHNOLOGY

SR515 (Wavelength: 515-950 nm)  
SR580 (Wavelength: 580-950 nm)



## PHOTOREJUVENATION

IPL technology delivers direct bursts of energy onto targeted areas of the skin to effectively treat pigmented and vascular lesions.

HR650 / HR650 XL (Wavelength: 650-950 nm)  
HR690 / HR690 XL (Wavelength: 690-950 nm)



## HAIR REMOVAL

IPL technology destroys hair follicles by targeting the pigment. XL applicators with larger spot sizes allow for faster and more effective treatments.


AC Dual (Wavelength: 415-480 nm)  
(Wavelength: 630-950 nm)



## ACNE TREATMENT

The AC Dual IPL applicator combines blue and red light delivered simultaneously with each pulse. The blue light targets and destroys acne-causing *P. acnes* bacteria, while red light reduces inflammation and promotes faster healing.


ALL THESE  
IPL APPLICATORS INCLUDE:  
**SMARTPULSE™  
TECHNOLOGY**



SmartPulse™ ensures precise and consistent energy delivery to the chosen target. It provides advanced pulse optimization for maintaining a perfect energy and wavelength profile throughout each pulse.

## REAL-TIME COOLING SYSTEM

Our innovative patent-pending cooling system monitors IPL applicator temperature 1,000 times a second for unparalleled safety and comfort. No downtime for cooling improves your productivity.



# RF APPLICATORS & OUR TECHNOLOGY



## NANOFRACTIONAL™ RADIO FREQUENCY

State-of-the-art patented tip technology with up to 700 pulses. The applicator also provides a therapeutic depth of penetration (up to 500 microns). The NanoFractional™ RF applicator with SmartScan™ provides varying energy density to enable ablation of the epidermis and coagulation of the dermis area, resulting in skin resurfacing with minimal discomfort.



## SMARTSCAN™ TECHNOLOGY

NanoFractional™ RF with SmartScan™ technology is delivered through a smaller footprint per pin (150 x 20 microns). The applicator features 160 pins per tip with 62mj per pin. SmartScan™ uses a unique algorithm to generate a customized pattern that enables maximum flexibility and control during the treatment.



# (MP)<sup>2</sup> APPLICATORS & OUR TECHNOLOGY

## DIAMONDPOLAR™ AND OCTIPOLAR™

The DiamondPolar™ applicator is used for the face and neck, and the OctiPolar™ applicator is used for the body. Both are powered by our patented (MP)<sup>2</sup> technology, which combines Multi-Polar Radio Frequency and Pulsed Electro-Magnetic Fields.



## MULTI-POLAR RADIO FREQUENCY

The Multi-Polar Radio Frequency delivers faster, more homogeneous heating to multiple tissue depths resulting in superior clinical outcomes with no pain or downtime for the patient.

## PULSED ELECTRO MAGNETIC FIELDS

Delivered simultaneously with Multi-Polar Radio Frequency, Pulsed Electro Magnetic Fields (PEMF) help to create more fibroblasts in the skin, which in turn boosts collagen. PEMF also promotes angiogenesis, thereby creating more pathways for blood flow and improving circulation to the skin.



# DELIVERING

# THE RESULTS



*Courtesy of Gitit Zucker, MD*  
Results after 7 treatments





*Courtesy of Venus Concept*  
Results after 2 treatments



*Courtesy of Megha Joy Shah, MD*  
Results after 5 treatments



*Courtesy of Gregory Antoniak, MD*  
Results after 2 treatments



*Courtesy of Tal Nachlieli, MD*  
Results after 1 treatment




*Courtesy of Rosenberg Plastic Surgery*  
Results after 6 treatments

# INNOVATION. PARTNERSHIP. SUCCESS.

Innovative Technology Powered by a Unique Business Model

When you choose Venus Versa™, you partner with Venus Concept. Our industry-unique business model ensures a higher return on your investment and further increases the success of your business.



# DEVICE SPECIFICATIONS

## TECHNICAL SPECIFICATIONS

Input Voltage	100-240 VAC , 12.5A, 50-60Hz, Single Phase
Dimensions	65 cm x 55 cm x 110 cm / 25.6 in x 21.7 in x 43.3 in (DxWxH)
Weight	35 Kg / 77 lb.

## SR515 / SR580 IPL APPLICATORS

Wavelength	515-950 nm, 580-950 nm
Optical Energy	5-25J/cm <sup>2</sup>
Pulse Duration	5-20 ms
Spot Size	10 mm x 30 mm
Pulse Repetition Rate	Up to 3 Hz

## HR650 (XL) / HR690 (XL) IPL APPLICATORS

Wavelength	650-950 nm, 690-950 nm
Optical Energy	5-20J/cm <sup>2</sup>
Pulse Duration	20-50 ms
Spot Size	10 mm x 30 mm (XL: 20 mm x 30 mm)
Pulse Repetition Rate	Up to 3 Hz

## AC DUAL IPL APPLICATOR

Wavelength	415-480 nm, 630-950 nm
Optical Energy/Fluence	5-25J/cm <sup>2</sup>
Pulse Duration	10-15 ms
Spot Size	10 mm x 30 mm
Pulse Repetition Rate	Up to 3 Hz

## NANOFRACTIONAL™ RADIO FREQUENCY APPLICATOR

RF Frequency	460 KHz
Fuse	3.15A, 250V
Max. Output Energy	62mJ/pin
Depth Ablation	500 microns
Number of Pins	160

## DIAMONDPOLAR™ & OCTIPOLAR™ APPLICATORS

RF Frequency	1 Mhz
RF Output Power	0-150 Watts
Number of (MP) <sup>2</sup> Synthesizers	Up to 8
Magnetic Pulse Frequency	15 Hz
Pulsed Magnetic Field Flux	15 Gauss

**Corporate Office:** 255 Consumers Road, Suite 110 | Toronto, Ontario, M2J 1R4

**USA Office:** 1880 N Commerce Pkwy., Ste 2 | Weston, FL 33326

**Tel: 888-907-0115 Visit us at [venusconcept.com](http://venusconcept.com)**

